

# DATASHEET

## Fiber Collimators



- SMA905 or FC connectors.
- Focal length of the collimator lens is  $F = 12\text{mm}$ .
- Operating wavelength is 200nm-1500nm.
- Can be used independently for collimation of fiber luminescence or in combination with various sample holders for spectral measurement of samples.

## Applications



Spectral analysis



Material inspection



Online inspection



R & D in science and technology



Customization etc.

All Rights Reserved: Glydden Science Pte Ltd. We have the rights to change designs and specifications of products without prior notice.